**PRIDE**<sup>®</sup>  
MOBILITY

**Jazzy**  
ASSIST

CE UK CA Copyright © 2026  
CMG001/Rev B/April 2026

**Owner's Manual**

# Contents

---

<b>Introduction</b> .....	<b>01</b>
<b>Safety Precautions</b> .....	<b>01</b>
1. Symbol Explanation	
2. Intended Purpose	
3. Intended Population	
4. Electrical Appliance Safety Precautions	
5. Product Operating Conditions	
6. Instructions Before Use	
7. Precautions for Use	
8. Excerpts from EMC-related Guidelines and Statements	
<b>Product Introduction</b> .....	<b>06</b>
1. Main Structure of the Product	
2. Main Technical Specifications	
3. Storage Conditions	
4. Safety Signs on the wheelchair	
<b>Use of wheelchair</b> .....	<b>08</b>
1. Unfolding and Folding of the wheelchair	
2. Introduction to the Controller	
3. Left and Right Armrest Controller Interchange	
4. Battery Disassembly and Assembly	
5. Switch to Push Mode	
6. Charging Introduction	
7. Connect to the App	
8. EZ-Pass Function Introduction	
9. Reversing Camera Function	
<b>Maintenance and Upkeep</b> .....	<b>16</b>
1. Service life	
2. Disposal	
<b>Cyber security</b> .....	<b>17</b>
<b>Wireless technology and related settings</b> .....	<b>17</b>
1. Wireless functions description	
2. Wireless security measures	
3. Information addressing wireless issues	
4. Information about any wireless coexistence issues and mitigations	
5. RF wireless communications information	
<b>Troubleshooting</b> .....	<b>18</b>
<b>Serious incident reporting</b> .....	<b>20</b>
<b>Packing List</b> .....	<b>20</b>
<b>Warranty</b> .....	<b>21</b>
<b>Symbols used on labels and packaging</b> .....	<b>22</b>

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---

Dear User:

Thank you for choosing the Jazzy Assist wheelchair (hereinafter referred to as the “wheelchair”). This manual is an important document that provides product information, operating instructions, and after-sales service information. Please read this manual carefully before use to ensure correct operation and maintenance. This manual contains information protected by patents and copyrights and may not be reproduced in whole or in part without the express written permission of our company. This manual provides operating instructions, limitations, warnings, and maintenance information in accordance with ISO 7176-15. Please read, follow, and retain this manual for future reference. If the product configuration, accessories, or software updates differ from this manual, please contact an authorized dealer for the latest information. The latest version shall prevail.

This wheelchair is a medical device. For the European Union (EU/EEA), it complies with Regulation (EU) 2017/745 (MDR) and bears the CE mark accordingly. For the United Kingdom, it complies with the UK Medical Devices Regulations 2002 (SI 2002 No. 618, as amended) and bears the UKCA mark accordingly.

The contents of this manual are based on the expectation that a mobility device expert has properly fitted the power chair to the user and has assisted the prescribing healthcare professional and/or the authorised Pride Dealer in the instruction process for the use of the product.

## Safety Precautions

---

### 1. Symbol Explanation

#### **Warning**

Text marked with this symbol indicates a hazardous situation that, if not avoided, may result in serious injury.

#### **Caution**

Text marked with this symbol indicates a hazardous situation that, if not avoided, may result in minor or moderate injury and property damage.

### 2. Intended Purpose

The intended use of this wheelchair is to provide indoor and outdoor mobility for people with mobility impairments who are able to operate a powered wheelchair safely.

### 3. Intended Population

This product is intended for people with mobility impairments who have difficulty walking long distances and who are able to operate a powered wheelchair safely.

## 4. Electrical Appliance Safety Precautions

### Caution

This wheelchair has been tested in accordance with ISO 7176-9 climatic test requirements.

- This test ensures that the user or caregiver has sufficient time to move the wheelchair out of the rain and into a safe location;
- Do not operate the wheelchair during thunderstorms;
- Do not expose this wheelchair to heavy rain;
- Do not leave this wheelchair in a damp place for an extended period of time;
- Do not use automatic car washes to spray or clean the wheelchair. If directly exposed to rain, smoke, or humid air, the wheelchair may experience electrical or mechanical failures, or rust.

## 5. Product Operating Conditions

### Caution

- Temperature: 0 - 40°C
- Relative humidity: 10-90%
- Atmospheric pressure: 78.9~106 kPa;
- No strong electromagnetic interference

## 6. Instructions Before Use

- Do not operate the wheelchair until you have carefully read and understood this user manual.
- Please refer to the technical specifications for the maximum load capacity. Overloading is strictly prohibited. This product is for single-person use only; do not carry passengers.
- Do not operate the wheelchair when you are under the influence of alcohol or when fatigued.
- Do not drive at night under low visibility conditions.
- The wheelchair must not be operated independently by users with visual impairment.
- Before use, please check whether the product has any abnormalities, including the battery, tires, joystick, and other components.
- This product has a certain level of waterproof ability (IPX4 rating), but it must not be stored or used in rainy or snowy weather.
- Turn off the wheelchair before getting on or off.

## 7. Precautions for Use

### Warning

Please carefully read the following precautions prior to use. Residual risks remain despite the implementation of risk control measures and verification / validation activities. Users must be aware that the wheelchair, like all mobility aids, **cannot eliminate all risks associated with powered mobility and use in real-world environments**. Safe use depends on correct operation, appropriate user capability, and suitable environmental conditions.

- The wheelchair is suitable for flat and level ground. Do not use it for extended periods on muddy, wet, uneven, or other unsuitable terrains.

- Avoid driving on roads with poor conditions or dangerous roads (such as grasslands, gravel roads, uneven surfaces, or bridge decks without guardrails). Always operate with caution and prioritize safety. Slow down when going uphill or downhill.
- Try to avoid going to crowded areas with dense pedestrian or vehicle traffic as much as possible.
- Always wear the seat belt while driving, especially on slopes, slippery surfaces, around obstacles, when turning, and in crowded environments.
- Before transferring into the wheelchair, raise the footrest. Do not stand on the footrest while transferring, as this may cause instability or injury.
- This wheelchair can be used on slopes up to 7°. Do not attempt to turn, accelerate sharply or brake sharply on slopes. Do not park on, or transfer on or off, slopes steeper than 23.8°. Slow down when turning on slopes to avoid rollovers. Always go up and down slopes in a straight line, not diagonally. Braking distance on slopes may be much longer than on flat ground, so please slow down when going downhill. Even when travelling on surfaces with a slope of <7°, there is still a risk of tipping over if the device is operated improperly.
- This wheelchair can cross steps or thresholds less than 2.5 cm in height. For steps or thresholds higher than 2.5 cm, direct crossing is prohibited. Please detour or use auxiliary facilities. Depending on the surface condition, approach angle, speed, load, and operating method, the wheelchair may not be able to cross a 2.5 cm step or threshold in a single attempt.
- When the wheelchair is in manual push mode, do not park it on a slope. Do not sit on the wheelchair at this time, to prevent accidents caused by it rolling downhill.
- Do not switch from drive mode to manual mode while the wheelchair is moving.
- Before use, be sure to check the battery level to prevent the wheelchair from stopping due to insufficient power.
- Once signs of instability (bumping, tilting, skidding) or any abnormal conditions (alarms, noises, vibrations, steering instability) occur, you should immediately slow down, pull over, turn off the power, and conduct a safety check.
- Do not operate the wheelchair outside its rated range (e.g., motor vehicle lanes, deep gravel pits, soft soil, icy roads, steep slopes, stairs, escalators, etc.).
- When operating a wheelchair, please place your feet firmly on the footrests to avoid the risk of your feet being caught in the wheels. Exercise extra caution when wearing loose or long clothing with a dangling belt, as these may get caught in the wheels.
- Do not connect any electrical or mechanical equipment to the wheelchair unless it has been officially approved. Using unauthorized accessories may void the warranty and pose safety risks.
- When lifting or transporting the wheelchair, be mindful of its weight and do not lift the device by any moving parts.
- Do not modify the wheelchair in any way that could compromise its safety.
- Do not drive the wheelchair in the middle of the road or in traffic lanes intended for motor vehicles. Be careful when crossing the road and watch out for pedestrians and vehicles.

- When transferring to a small vehicle for travel, it is strictly prohibited for the user to sit on the wheelchair. Please transfer the user to the vehicle's original seat and use the vehicle's seat belt.
- Do not bring the wheelchair into the MRI examination room to avoid causing personal injury or equipment damage.
- Do not make a U-turn in narrow spaces near edges (steps, potholes, etc.).
- Always ensure that the anti-tilt wheels are in good condition. Mounting heavy objects on the wheelchair will cause the center of gravity to shift, so please slow down and operate with extra caution.
- Use low-speed mode and accelerate slowly in indoor, crowded, or complex environments. Prolonged high-speed driving is prohibited.
- Please avoid high-speed collisions between the wheelchair and obstacles or walls. It is prohibited to use the wheelchair for climbing or to bear concentrated point loads (such as standing only on the footrest or hand-rail). When transporting or crossing slopes/doorsteps, caregivers must provide assistance from a stable position to prevent structural damage caused by torsional stress.
- Long-term exposure to sunlight, proximity to heat sources, or continuous high-current operation may lead to overheating. If the temperature of the housing or motor is too high, use should be stopped immediately, and inspection should be carried out after it has cooled down. Low temperatures will reduce battery and braking performance. Before starting or charging in cold conditions, the wheelchair should be placed in a dry and ventilated environment to return to room temperature. Any abnormal conditions (burning odor, abnormal noise, sudden performance drop) require immediate power-off and contacting the service provider.
- Regularly clean and disinfect components that come into contact with the human body (seat cushion, armrests, joystick) as necessary, following the cleaning methods in this manual. Keep the wheelchair dry and well-ventilated. Do not use high-pressure water, steam cleaning, or soaking. Avoid liquid from entering electronic components or connectors. If contaminated by body fluids, immediately follow medical cleaning procedures and wear appropriate personal protective equipment.
- Do not allow others to pull or drag the wheelchair, and do not use it as a towing device.
- Do not leave the wheelchair—especially the battery—below 0°C for extended periods (e.g., overnight). Store it indoors in a dry, warm place. If stored in an unheated garage or other cold area, remove the battery and keep it indoors.
- Never charge in cold conditions. Move the wheelchair or battery indoors and wait until the battery warms to above 5°C (ideally near room temperature) before charging.
- When checking in the wheelchair for air travel, use the approved wheelchair transport bag. Before check-in, fold the front caster wheels. For air transport, large aircraft are recommended as a priority.
- Please do not press on the armrest screen area when boarding or alighting, to avoid unnecessary damage.

## 8. Excerpts from EMC-related Guidelines and Statements

This product complies with **IEC 60601-1-2:2014/AMD1:2020** standard, meeting the requirements in both electromagnetic immunity and emissions.

External electromagnetic interference (such as radio, television transmitters, CB radio waves, garage door openers, cordless phones, etc.) may cause the following problems:

- Brake failure, automatic startup, or abnormal steering;
- The control system is permanently damaged;
- Unexpected movement, performance degradation, or loss of control.

If the above anomalies occur near a strong electromagnetic field source, immediately release the joystick/turn off the power; move the wheelchair to a safe distance before attempting to use it again.

Keep a safe distance away from the following strong electromagnetic radiation sources:

- Walkie-talkies, industrial transmitters, substations;
- MRI (Magnetic Resonance Imaging), RF treatment equipment;
- Radio Frequency Identification (RFID) readers, metal detectors, electronic article surveillance systems;
- Near Field Communication (NFC) systems, Wireless Power Transfer (WPT);
- Electrocoagulators, electrosurgical devices, diathermy devices.


### **Warning**

- Modifying the product or using electrical accessories not approved by the manufacturer is strictly prohibited, as it may increase electromagnetic radiation or reduce anti-interference ability.
- Avoid using this device adjacent to or stacked with other devices; if it is necessary to do so, observe whether both are operating normally.
- Only use accessories, sensors, and cables specified or provided by the manufacturer; non-specified accessories may degrade EMC performance and cause abnormal operation.
- If the wheelchair is used within 1.5 km of an AM/FM/TV broadcast antenna, verify that it operates normally before use.
- When the AC input voltage is interrupted, the device will stop charging the battery; after the power is restored, charging will automatically resume (this performance degradation is acceptable, with no safety risks or loss of basic performance).

The above content must be strictly adhered to to ensure the safe and stable operation of the equipment.

# Product Introduction

## 1. Main Structure of the Product

## 2. Main Technical Specifications

Item	Specification
Net weight (excluding battery)	15.9 kg (35.05 lb)
Maximum load	136 kg (300 lb)
Overall dimensions (L×W×H)	1035 × 620 × 872 mm (40.75 × 24.41 × 34.33 in)

Folded dimensions (L×W×H)	820 × 620 × 380 mm (32.28 × 24.41 × 14.96 in)
Seat width	460 mm (18.11 in)
Seat depth	437 mm (17.20 in)
Backrest height	438 mm (17.24 in)
Maximum speed	6 km/h (3.73 mph)
Obstacle clearance height	≤ 25 mm (0.98 in)
Level-surface braking distance	≤ 930 mm (36.61 in)
Climbing ability	≤ 7°
Minimum turning radius	852 mm (33.54 in)
Theoretical travel distance	10Ah: with 100 kg load, up to 13.6 km (8.45 mi) 20Ah: with 100 kg load, up to 27.2 km (16.90 mi)
Battery capacity	10Ah (standard) / 20Ah (optional)
Front/Rear wheel size	7" / 10"
Tire type	Solid rubber tire (front), pneumatic rubber tire (rear)

Note:

1. The weight-bearing capacity is tested using a dummy or a person of the same weight.
2. The maximum cruising range is tested on flat roads at an average speed, and the actual cruising range may vary due to factors such as road surface and driving habits.

### 3. Storage Conditions

- Packaging storage temperature: -20°C to 65°C;
- Relative humidity: ≤ 90%;
- Atmospheric pressure: 56–106 kPa;
- Store indoors in a well-ventilated area free of corrosive gases. Do not place heavy loads on the package;
- Maximum stacking height (packaged): 4 units;
- Do not store the product near heat sources or in direct sunlight for extended periods;
- If the product is stored for a long period without use, fully charge it, turn off the power switch, and then store it;
- If stored for more than one month, fully recharge the battery before continuing to store or use it.

## 4. Safety Signs on the wheelchair

The following symbols are used on the wheelchair to indicate warnings, mandatory actions, and prohibited operations. It is very important to read and fully understand these symbols.

### Pinch Point Warning:

Do not put your hands into moving mechanical parts to avoid injury.

### LiDAR Cover Warning:

Non-load-bearing component, do not step on it.

### Footrest Platform Warning

Do not stand on the footrest, as this may cause tipping or injury.


When the brake lever is set to the **[Ride]** icon, the wheelchair is in **electric-drive mode**. When it is set to the **[Push]** icon, the wheelchair is in **push mode**. When switching modes, place the wheelchair on a flat surface and stand to one side while operating the brake release switch.

## Use of wheelchair

### 1. Unfolding and Folding of the wheelchair

#### Unfold


- 1 Step on the folding footrest, lift the handle upwards and push it forward

- 2 Open the backrest

- 3 Lift the armrests on both sides


- 4 Lower the footrest

## Fold

- 1 Slightly lift the armrest, press and hold the lock release button, and lower the armrests on both sides


- 2 Step on the folding pedal, lift the backrest handle and then release your foot



- 3 Pull back the locking link and fold the backrest

- 4 Retract the footrest (optional)


- 5 Pulling the locking mechanisms on both sides of the front support wheel allows the front support wheel to be bent backward (not mandatory)

Caution: Please pay attention to the pinch point warning label on the wheelchair and avoid pinching your fingers.

## 2. Introduction to the Controller

## Directional Joystick

When the wheelchair is powered on, move the joystick to control speed and direction. Pull the joystick backward to display the rear-view camera image on the screen.

## Power/Lock

Short press to power on. Press and hold for 3 seconds to power off. When powered on, short press to enter lock mode; the joystick is disabled.

## EZ-Pass

Short press to enable this mode. When enabled, it makes minor steering corrections to help wheelchair users pass through narrow doors or passages more easily. Short press again to disable this mode.

## Horn

Press to sound the horn.

## Light

Short press to turn the lights on or off.

## SOS

Press and hold for 3 seconds to sound the alarm. Press and hold for 5 seconds to send an SMS message to the emergency contact (Only effective when the wheelchair is connected to the mobile app).

## Menu

Short press to enter the menu. Use the [speed up]/[speed down] buttons to switch pages. In the menu, press and hold the [speed down] button to display the exit prompt, then short press the [Menu] button to exit.

## Speed Up

Short press to increase by one level, up to a maximum of Level 5. In the menu, short press to go to the next page. On the volume adjustment screen, short press to increase volume.

## Speed Down

Short press to decrease by one level, down to Level 1. In the menu, short press to go to the previous page. On the volume adjustment screen, short press to decrease volume.

## Push Mode

Press and hold the [Push] button for 3 seconds to release the brakes so the wheelchair can be pushed manually. Short press the [Push] button or move the joystick to return to electric drive mode. You can change the button function in the app.


### 3. Left and Right Armrest Controller Interchange

1. Remove the battery before swapping the armrests.
2. Use a hex wrench to remove the four hex socket screws under the armrest controller, then disconnect the wiring harness connecting the left and right controllers.
3. Swap the left and right controllers and reconnect the wiring harness. Insert the controller back into the slot and tighten the hex socket screws to complete the replacement.

### 4. Battery Disassembly and Assembly

#### Remove the battery

- 1 Lower the battery handle.

- 2 Hold the handle and pull the battery out.

#### Install the battery


- 1 Secure the soft rubber plug at the battery charging port.


- 2 Lift the battery handle.

- 3 Slide the battery all the way along the chute and ensure it is fully seated.

- 4 Reset the battery handle.

## **Caution**

- When turning on the wheelchair for the first time after unpacking, press and hold the button on the battery for 3 seconds to activate it. A green light indicates successful activation. You can then turn it on and use the wheelchair.

- If you cannot push the battery all the way in, check whether the soft rubber plug on the battery charging port is secured properly.

## **5. Switch to Push Mode**

Method 1: Use the brake release switches. First, make sure the wheelchair is locked or powered off, then pull down both brake release levers to enter push mode.

Method 2: Press and hold the [Push Mode] button for 3 seconds to switch to push mode, and the wheelchair can be pushed manually. Move the joystick or short press the [Push Mode] button to return to electric drive mode. The wheelchair will move only when the joystick is moved.


## **Caution**

- When the wheelchair is powered on and not in lock mode, toggling the brake release lever will cause a fault state, and the joystick will be unable to control movement.
- Both brake release switches must be in the driving position. Otherwise, the wheelchair will enter fault mode.

## 6. Charging Introduction

### Charging Method

- Method 1: Connect the charger. Align the red dot on the plug with the red dot on the charging port under the armrest, and insert it in the same direction to start charging.

- Method 2: Remove the battery. Align the red dot on the plug with the red dot on the battery charging port, and insert it in the same direction to start charging.

### Charging Status Description

1. When charging while the device is powered on, the screen will display the charging level in real time.
2. When charging while powered off, the battery indicator will flash alternately to indicate charging. The indicator will turn off when the battery is fully charged.

#### **Caution**

- The device cannot be used while the battery is charging.
- When the charger is powered on and connected to the wheelchair for more than 12 hours, it will automatically cut off power and enter overcharge protection. Disconnect power to exit overcharge protection; only then can it resume charging the wheelchair.
- Store the battery in an environment above 0°C. Before storage, make sure the battery is fully charged. Recharge the battery at least once per month.

- The battery may become over-discharged and unable to accept a charge.
- Foreign objects may cause a short circuit.
- Read this manual before using the battery charger.
- The battery charger is for indoor use only and must not be exposed to rain or other wet conditions.
- Explosive gases may be generated during charging, so the wheelchair and battery charger should be kept away from sources of ignition, such as flames and sparks.
- Use only the official charger to charge the battery.

## 7. Connect to the App

1. Scan the QR code to download and install the companion app for the wheelchair.

2. Open the app, then follow the in-app instructions to add the device.
3. For more information on additional features (such as SOS emergency assistance and app remote control), please refer to the app.

### Caution

- The companion app for this product supports Android 14 and iOS 16.6 or later.
- This product supports only 2.4 GHz Wi-Fi networks.
- Due to updates, actual operations may differ slightly from the above description. Please follow the in-app instructions.
- This product can be bound to only one mobile app account at a time. You can bind it to another account only after unbinding it first.

## 8. EZ-Pass Function Introduction

### EZ-Pass Unlock Method

### Caution

This feature is off by default, and the EZ-Pass button on the wheelchair is disabled. You must enable it after the wheelchair is bound to the mobile app. By default, this feature is limited to Speed Level 2. The speed limit can only be removed after completing the in-app safety confirmation steps.

When EZ-Pass is activated, the wheelchair will proactively make minor steering adjustments to provide steering assistance in narrow passages. When it determines that it cannot pass safely, it will slow down and stop to reduce the risk of collision or tipping.

This feature can reduce the probability of collisions, but it cannot completely avoid them. Do not rely solely on this feature, especially in certain special situations. Users can choose whether to enable this feature; it is only an auxiliary passage function.

## Objects and conditions that EZ-Pass may not detect

Category	System Capability Boundary	Typical Scenario
Height Direction Target	Unable to detect low objects with a height $\leq 26.4$ cm	Items on the floor with a height less than 26.4 cm, such as children's toys (toy cars/strollers), slippers, floor mats, etc.
	Unable to detect suspended objects with a clearance above the ground $> 26.4$ cm	Ceiling-mounted shelves or racks, decorative lighting fixtures, chandeliers, etc., with their lowest edges higher than 26.4 cm
Small/Fine Target	Unable to detect objects with a diameter or width $< 2$ mm	Wires, thin rods, narrow cables
Transparent Target	Unable to detect transparent objects (objects with visible light transmittance $> 85\%$ and near-infrared transmittance $\geq 90\%$ )	Glass doors, acrylic sheets, transparent plastic films
High reflectivity target	Unable to detect high reflectivity objects (reflectivity $> 85\%$ )	Mirrors, stainless steel components, polished metal panels
Low reflectivity/ absorptive target	Unable to detect black or highly absorbent objects (near-infrared reflectance $< 5\%$ )	Black rubber, matte black ABS plastic, dark carbon fiber, black fabric, coal
Road surface depression/level difference	surface depressions, drops, or slopes with a height difference $< 26.4$ cm	Sunken road, downhill slope, descending steps

### Warning

Even if an obstacle falls within the detectable categories, the LiDAR may still fail to detect it. Please make sure you understand how this function works, including where it can assist and where it may not detect obstacles reliably. If you have any questions, please do not activate this function to ensure your safety. Any risks arising from use beyond the intended purpose (including misuse, collision, or injury) are not within the manufacturer's liability.

### Caution

It is recommended that users temporarily turn off the EZ-Pass function when the following situations occur:

- Rain/Fog/Condensation: Scattering/attenuation caused by water droplets or fog droplets in the air;
- Strong direct infrared radiation: The sensor window is directly exposed to sunlight or a strong infrared source;
- Severe stains/scratches: contamination, scratches on the LiDAR cover surface, wrinkling of the film, etc.;
- Scratches, dust, or dirt on the LiDAR cover can reduce detection performance and should be cleaned promptly;
- When passing through doors that require manual pushing, such as rolling doors and revolving doors.

## 9. Reversing Camera Function

### Activating the reversing camera

Pull the joystick backward. The screen will display the rear view and provide alerts for detected pedestrians, while the wheelchair moves backward at the same time. After stopping, the image will remain on the screen for 3 seconds and then disappear. Pushing the joystick forward will clear the image immediately.

### Hazard detection warning

When reversing, the rear camera may detect pedestrians and highlight them on the screen. If a pedestrian is detected within 1.5 m, an audible warning will sound.

## Maintenance and Upkeep

---

- Regularly clean dust from sensor surfaces (e.g., LiDAR sensors and the rear camera) to ensure proper operation of all functions.
- With power off, check the joystick daily for damage and ensure it returns correctly to the neutral position. If any defects are found, contact the distributor immediately. Do not attempt to repair it yourself.
- Check the seat, armrests, and footplate daily for loose parts, excessive wear, or damage. If any issues are found, contact the distributor.
- Visually inspect tire treads regularly. Replace worn tires at an authorized distributor if necessary.
- Check rear tire pressure periodically; recommended pressure is 35 PSI / 2.4 Bar. Inflate if below this value. Low pressure accelerates tire aging and affects range.
- Inspect all exposed wiring harnesses and connectors for wear, damage, or corrosion. Any damaged wires or connectors must be repaired or replaced by the distributor.
- Keep all electronic components (e.g., charger and joystick controller) away from moisture. If they become damp, do not use them again until completely dry.
- Check all fasteners regularly for loosening, wear, or corrosion. If any abnormalities are found, turn off power immediately and contact the distributor.
- Avoid exposing the joystick controller to harsh weather. Moisture may damage the controller and void the wheelchair's warranty.

## 1. Service life

The expected service life of this product is 5 years from the date of manufacture. Do not continue to use the product beyond its intended service life.

## 2. Disposal

When disposing of the device, its parts and the battery or its parts, be sure to follow the local disposal guidelines. If there is anything that is unclear, consult your dealer or customer support.

# Cybersecurity

---

This wheelchair incorporates embedded control software and wireless communication functions.

To ensure safe operation:

Do not install or modify any software without manufacturer authorization.

Wireless data are encrypted (AES128). Pairing must be done safely.

In case of software error, the system enters safe stop mode.

Personal data (if collected) are protected and deleted upon request in accordance with GDPR.

# Wireless technology and related settings

---

## 1. Wireless functions description

The wheelchair can be controlled by mobile phone via Bluetooth 5.4.

## 2. Wireless security measures

At the application layer, interactive data is protected by hardware and cannot be read directly. If there are multiple instances of data errors or illegal content verification, the connection will be interrupted and the [offending party] will be added to a blacklist until the coverage scope allows for its removal. This ensures that signals are transmitted as intended.

## 3. Information addressing wireless issues

If the mobile phone cannot control the wheelchair correctly, stop using the remote-control function immediately and use the on-board controls instead.

Please remove other wireless devices from the environment, and connect mobile phones to wheelchair again, try controlling wheelchair with mobile phones. Do not use remote control function until the mobile phone can control the wheelchair correctly.

## 4. Information about any wireless coexistence issues and mitigations

When using the mobile phone, try to keep other wireless devices away from the wheelchair, at least 2 m away.

## 5. RF wireless communications information

Please take attention that changes or modification not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Operation is subject to the following two conditions:

(1) This device may not cause harmful interference.

(2) This device must accept any interference received, including interference that may cause undesired operation.

## Troubleshooting

When an abnormal prompt appears or the device cannot be used normally, please refer to the table below for an initial check. If the issue persists after following the suggested actions, stop using the device and contact your authorised dealer or customer support. This manual lists only basic troubleshooting methods for common abnormalities. For more detailed fault information, real-time status details, or further guidance, open the official app and go to **Device Diagnostics / Help & Support**.

Common Fault / Symptom	Possible Cause(s)	User Troubleshooting (in order)
Unable to power on	Battery depleted; Battery not activated; Battery not fully seated	① Charge for 30 minutes and try again.② Press the battery button briefly to check whether the indicator lights turn on. If not, press and hold for 3 seconds to activate the battery, then try again.③ Charge for another 30 minutes and try again.④ Check whether the battery is fully seated.⑤ If there is still no response, contact after-sales service.
Can move, but speed is significantly reduced	Low battery; Low temperature reduces battery output; Low tire pressure / high rolling resistance; Motor/controller overheating with power derating; Excessive load	① Fully charge the battery.② Move indoors and allow it to warm up before use in low-temperature conditions.③ Check tire pressure / check whether the wheels are jammed by debris.④ Stop and let the device cool for 10 minutes.⑤ Reduce load / climbing intensity.
Sudden shutdown during driving	Undervoltage protection (low battery); Overcurrent/overtemperature protection; Excessive road resistance	① Check the battery level. If low, charge immediately or replace the battery.② Let it cool down, then restart.③ Check whether any wheel is jammed.
Charger is plugged in but not charging	No power at the outlet; Charger malfunction; Poor contact at the charging port; Battery temperature too low (<5°C); Battery aging	① Check whether the battery has been activated.② Check the outlet / confirm the outlet has power.③ Unplug and replug the charging connector and the battery terminal.④ Ensure the battery temperature is above 5°C before charging.⑤ If still not charging, contact after-sales service to check the charger/battery.
Driving range drops sharply	Battery aging; Low temperature; Low tire pressure; Frequent hill climbing / grass or other high-resistance surfaces; Brake drag	① Warm up indoors and then re-evaluate the range.② After fully charging, perform a short-distance range test.③ Check whether the tires/casters roll smoothly and whether tire pressure is normal.④ Reduce slope / load.⑤ Replace the battery if necessary.
Vibration / pulling to one side / unstable steering	Uneven tire pressure; Loose wheel bearings; Tire wear; Uneven ground; Front-wheel shimmy	① Check that left and right tire pressures are the same.② Inspect tires for bulges or severe wear.③ Remove any entangled debris from the wheels.④ Test drive at low speed.⑤ If the issue remains obvious, contact after-sales service for alignment/bearing inspection.

Stops unexpectedly during driving	LiDAR sensor blocked by clothing; Large items mounted/hanging on the side	① Check whether the LiDAR sensor is blocked by clothing, blankets, etc.② Check whether side bags/accessories are blocking the LiDAR sensor.
Armrest lock button cannot be pressed during folding	Armrest is under downward pressure; Button malfunction	① Lift the armrest first, then press and hold the lock button and check whether the armrest can be lowered.② If still cannot be pressed, contact after-sales service.
Foldable footrest cannot be pressed down during unfolding	Not operated according to the manual; Foldable footrest malfunction	① While lifting the backrest/armrest, press down on the foldable footrest.② If still cannot be pressed down, contact after-sales service.
Charger/battery becomes abnormally hot or emits an odor	Poor contact; Charger malfunction; Internal battery abnormality; Poor ventilation	① Stop charging immediately and disconnect the power.② Move to a well-ventilated area and observe.③ Check whether the plug/connector is blackened or loose.④ Do not continue using; contact after-sales service.

## Serious incident reporting

---

In case of a serious incident (death, temporary or permanent serious deterioration of health, or serious public health threat) related to this device, Please contact the manufacturer or the EU/UK Authorised Representative. Keep records of incident details, device condition, and patient/user status. Failure to report serious incidents may compromise patient safety and regulatory compliance.

## Packing List

---

Item Name	Qty.	Remarks
Wheelchair	1	Includes battery
User Manual	1	—
Quick Start Guide	1	—
Safety Belt Kit	1	Includes safety belt, screws, and L-shaped Allen key
Charger	1	—
Cushion	1	—
Under-seat Mesh Basket	1	—

# Warranty

---

## LIMITED TWO-YEAR WARRANTY

Load-bearing frame components, including: platform, fork, seat posts, and frame welds.

The drive system, including the differential gearbox, motor, and brakes.

We provide a warranty of two (2) years from the date of purchase for all damage resulting from manufacturing and material defects. Your statutory rights remain unaffected by this warranty.

All electronic parts, including the controls, are covered by a two (2) year warranty. The controller and battery charger must be serviced by your authorized Pride/Quantum Rehab dealer. Any attempt to open or disassemble these parts will void the warranty.

## EXCLUSION OF WARRANTY

Excluded from the warranty are wear parts and parts that must be replaced due to wear and tear (tires, belts, lamps, upholstery, plastic covers, carbon brushes, fuses, and batteries), as well as product damage caused by improper use or accidents, for which Pride and its representatives assume no responsibility. This warranty does not cover labor costs or customer service.

## BATTERIES

The batteries are covered by a manufacturer's warranty of six (6) months.

Not covered by the warranty is any gradual decline in battery performance due to a continued discharged state, storage at low temperatures for extended periods, or wear resulting from heavy use.

## WARRANTY DISCLAIMER

Warranty service can be provided by your authorized Pride/Quantum Rehab dealer. For information about service costs, please contact your authorized Pride/Quantum Rehab dealer.

## REPLACEMENT UNITS

The availability of replacements is at the discretion of the supplier and not the manufacturer. For further information regarding replacement units, please contact your authorized Pride/Quantum Rehab dealer.

# Symbols used on labels and packaging

Symbol	Introductions	Symbol	Introductions
	Manufacturer Name Address		Refer to instruction manual
	Warnings and Precautions		CE marking and UKCA marking
	Medical device		The device's protection against splashing water from any direction without functional impairment.
	Keep dry		Keep away from Sunlight
	WEEE (Waste Electrical and Electronic Equipment) regulations require that waste products be disposed of in accordance with legal guidelines.		Unique Device Identification
	Maximum number of identical transport packages/items which may be stacked on the bottom package, where "n" is the limiting number.		This is the correct upright position of the distribution packages for transport and/or storage.
	BATTERY-POWERED VEHICLE		FRAGILE
	Dispose of Waste Properly		Recycle
	European Authorized Representative		UK Responsible Person

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Beijing OrionStar Technology Co., Ltd
SUNGO Europe B.V.
SUNGO Certification Company Limited
Pride Mobility Products Ltd.
Pride Mobility Products GmbH