

# Jazzy<sup>®</sup> EVO 613Li



**PRIDE<sup>®</sup>**  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039  
pridemobility.eu . info@prideeurope.nl

De Jazzy EVO 613 Li, uitgerust met lithiumbatterijen, belichaamt de filosofie van voortdurende verbetering en innovatie die wordt gedeeld door Pride Mobility®. Wij geloven in het verleggen van grenzen, zodat u nooit achterop raakt.



Compacte voetplaat



13" massieve banden



Lithium ion accu

## Specificaties

Maximale belasting	136 kg
Snelheid	Tot 8,2 km/u
Maximale actieradius	Tot 33,8 km
Draaicirkel	559 mm
Bodemvrijheid	63 mm
Lengte	889 - 1041 mm
Breedte	559 mm
Aandrijf wielen	13" massief
Zwenkwielen	5" massief
Gewicht zonder accu's	64,5 kg
Accu's	Type: Lithium ion
Meer?	Download het specificaties blad op de website

## Kenmerken

- Technologie voor middenwielaandrijving
- Actieve-Trac-vering
- 13" massieve aandrijfbanden
- Lithium ion accu
- Tot 8,2 km/u

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

## Beschikbare kleuren

